

# ERKA SWITCH 2.0 CONFORT ANEROID SPHYGMOMANOMETER - 1 tube



Product Code	32698
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1147713
CND	C9006
EAN/UPC	4250296150729
UMDNS	16-156
PARAF	971057195

## Description

ERKA SWITCH 2.0 CONFORT - dial Ø 56 mm - 1 tube grey - 293.204.93

It offers outstanding ergonomics and simple adaptation for left and right-handle users. It is manufactured from shock-resistant polycarbonate and has passed the shock test according to AAMI SP9. The sight glass is chemically stable and unbreakable. The brass valve screw works smoothly and is adjustable to exact settings.

MR compatibility. Measuring range 0-300 mmHg  $\pm$ 1%. Supplied with pouch. Made in Germany. 4 year warranty hand held, 2 year manometer.