



Product Code	32702
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1711340
CND	C9006
EAN/UPC	6970384251876
UMDNS	16156
PARAF	980925820

Description

SHOCK RESISTANT SPHYGMO - nylon-cuff

High quality, shock resistant sphygmomanometer matching the ISO 81060 standard.

One piece design with damped system manometer

Creative design with a new structure that fixes all parts together closely to make it more steadier. The damped system absorbs impacts.

Other features

- large dial Ø 55 mm easy to read
- high precision air-release valve, provides precise rate adjustment and quick deflation
- adult nylon cuff with latex-free bladder
- soft bulb for fast inflation.

Multilanguage box and manual: GB, FR, IT, ES, PT, DE.