



Product Code	32703
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	18158
CND	C9006
EAN/UPC	8023279327038
UMDNS	16156
PARAF	903174718
GTIN	08023279327038
GMDN	16156

Description

YTON ANEROID SPHYGMOMANOMETERS

YTON - with incorporated - stethoscope - black cuff

Economic sphygmomanometers with nylon cuff with velcro closure, reliable manometer with clip to be fixed on the cuff. Scale 0-300 mmHg. Air release at closed tap max 4 mmHg/minute. Manual setting of deflation possible up to ? mmHg/sec. From 260 mmHg to 15 mmHg max deflation time 10 sec. Scale graduated every 2 mmHg. Floating zero (the pointer has no stop point but swings freely), nylon velcro cuff, latex bulb with chromium plated valve. Chromium-plated air tap with air release regulation by screw valve. Supplied with a simulated-leather soft pouch with zip in a carton box.

Available in 5 versions: standard, with one piece cuff, latex free, with incorporated stethoscopes or kit (sphygmomanometer+stethoscope).

Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, Arabic.