

SIRIO PEDIATRIC KIT with 3 pediatric cuffs



| | |
|--------------|-------|
| Product Code | 32705 |
|--------------|-------|

| | |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

| | |
|---------------|---|
| Minimum order | 1 |
|---------------|---|

| | |
|------|------------------------|
| Type | Kit of medical devices |
|------|------------------------|

| | |
|-------|--|
| Class | |
|-------|--|

| | |
|--------|----|
| UK-REP | No |
|--------|----|

| | |
|--------|----|
| CH-REP | No |
|--------|----|

| | |
|---------|---------------|
| EAN/UPC | 8023279327052 |
|---------|---------------|

| | |
|-------|-----------|
| PARAF | 987470212 |
|-------|-----------|

Description

SIRIO PEDIATRIC KIT - sphygmo with 3 pediatric cuffs

Sirio with 3 pediatric cuffs: minicuff (32886) 4 months - 3 years, minicuff (32888) 3-7 years, minicuff (32880) 7-10 years in a nylon bag.