



Product Code	32709
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	18161
CND	C9006
EAN/UPC	8023279327090
UMDNS	16156
PARAF	903174744
GTIN	08023279327090
GMDN	16156

Description

ROMA SPHYGMO - pediatry - 1 tube

With pediatric cuff 35 x 11 cm.

Classic, lightweight, reliable palm model aneroid sphygmo with certified manometer ABS housing, integral PVC bulb and air release valve for one-hand use, provided with handle, blue nylon, velcro closure, calibrated cuff and 1 tube TPU bladder.

Easy to read: large dial Ø 60 mm, with rubber protective ring.

Supplied in a nice nylon carrying bag with external velcro pocket and cardboard box.

Latex free.