

BABYPHON PEDIATRIC SPHYGMOMANOMETER with 3 cuffs - metal - 1440



Product Code	32710
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2630625
CND	C9006
EAN/UPC	8023279327106

Description

BABYPHON® PEDIATRIC SPHYGMO with 3 cuffs - metal - 1440

Palm pediatric sphygmomanometer with 3 cuffs

- fast inflation of the cuff due to ergonomic spoon-shaped stainless steel handle
- easy reading metal manometer with aluminium scale Ø 49 mm
- lifelong precision due to almost non-aging precision movement with specially tempered copper-beryllium diaphragm for pressure-loading capacity up to 600 mmHg
- linear scale up to 300 mmHg, ± 3 mmHg
- bear imprint on the cuffs

Supplied with 3 pediatric cuffs in a padded nylon zipper case:

Newborn: size 15.5x5 cm - circumference 5-8 cm

Infant: size 23x7.5 cm - circumference 8-13 cm

Baby: size 35.5x10 cm - circumference 13-17 cm

Made in Germany.