



| | |
|---------------|----------------|
| Product Code | 32711 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I M |

| | |
|------------|---------------|
| UK-REP | Yes |
| CH-REP | Yes |
| RDM (NSIS) | 101769 |
| CND | C9006 |
| EAN/UPC | 4045396001817 |
| PARAF | 987475478 |

Description

MINIMUS® II sphygmo

Minimus® II are high quality, compact and corrosion resistant aneroid sphygmomanometers, featuring a chromium-plated screwed tube connector and ergonomic stainless-steel spoon-shaped grip.

Latex free bulb with chromium-plated air release valve.
Precision air release valve with fine adjustment and wearfree.
Linear aluminium graduated-scale Ø 49 mm, up to 300 mmhg.
Maximum error tolerance: ± 3 mmhg.

Supplied with adult cuff and vinyl zipper case.

Made in Germany.

Metal connectors for easy and fast changes of cuffs