

MINOR 2 SPHYGMOMANOMETER - nylon velcro cuff



Product Code	32714
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1860188
CND	C9006
EAN/UPC	8023279327144
UMDNS	16156
PARAF	922361124

Description

MINOR 2 SPHYGMO - 2 tubes - black nylon velcro cuff

High quality sphygmomanometer with traditional design and modern features.

- high precision manometer, adjustable by turning a screwdriver, for perfect calibration any time
- 2 tubes system for quick inflation
- Available with washable coated fabric velcro cuff or cotton hook cuff

Manual in: GB, FR, IT, ES, PT, DE, GR, Arabic

Adjustable gauge for perfect calibration

LATEX FREE