

R1 SHOCK PROOF SPHYGMOMANOMETER with velcro cuff - 1250-107



Product Code	32716
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396159846

Description

R1 SHOCK PROOF SPHYGMO with velcro cuff - 1250-107

Exact precision that lasts a lifetime.

- shock-proof up to a falling height of 120 cm
- valve base and air duct with metal precision valve
- two-component ABS casing
- metal threaded tube connection at the top of the casing for fast cuff exchange and ergonomical operation
- maximum error tolerance of ± 3 mmHg
- supplied in zipper bag

Non-resistant opening of air release valve ensures precise and comfortable operation

Scale: \varnothing 52 mm black with red fluorescent pointer

Fast, ergonomic inflation due to optimized bulb design with integrated spoon

Weight: 173 g

Made in Germany.