



Product Code	32718
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396002586
---------	---------------

Description

RI-SAN® SPHYGMO - green

- High quality palm-style sphygmomanometer with ergonomic design.
Air release valve conveniently regulated by finger-tip control.
- calibrated nylon cuff with velcro
 - latex-free bulb for pumping up the cuff
 - tube connection at the top of pressure gauge for unhindered blood pressure measurement
 - microfilter protects air release valve and measuring system
 - specially tempered copper-beryllium diaphragm almost non-aging
 - the precision movement guarantees reliable measurement
 - diaphragm can be pressurized up to 600 mmHg without decalibrating the sphygmo.

Easily readable scale Ø 64 mm, from 0 to 300 mmHg, ±3 mmHg
Precision air release valve with fine adjustment and wear-free
Suitable for right or left handed users