



Product Code	32722
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2934948
CND	C9006
EAN/UPC	4045396112902

## Description

RI-SAN® SPHYGMO with stethoscope - slate grey - 1444-141

High quality palm-style sphygmomanometer with ergonomic design.

Air release valve conveniently regulated by finger-tip control.

- calibrated nylon cuff with velcro
- latex-free bulb for pumping up the cuff
- tube connection at the top of pressure gauge for unhindered blood pressure measurement
- microfilter protects air release valve and measuring system
- specially tempered copper-beryllium diaphragm almost non-aging
- the precision movement guarantees reliable measurement
- diaphragm can be pressurized up to 600 mmHg without decalibrating the sphygmo.

Easily readable scale Ø 64 mm, from 0 to 300 mmHg, ±3 mmHg

Precision air release valve with fine adjustment and wear-free

Suitable for right or left handed users