



|               |                |
|---------------|----------------|
| Product Code  | 32725          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I M            |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 10049          |
| CND        | C9006          |
| EAN/UPC    | 8023279327250  |
| UMDNS      | 16156          |
| PARAF      | 903174922      |
| GTIN       | 08023279327250 |
| GMDN       | 16156          |

## Description

LONDON SPHYGMO - adult - black cuff

Professional, stylish, reliable, lightweight and latex free sphygmomanometer with ABS manometer, integral PVC bulb, high quality precision mechanism, velcro closure calibrated nylon cuff and one tube TPU bladder. Dial Ø 55 mm, with specially designed eye-catching colour printing and easy to read. Special outer ring protects the manometer from shocks and the glass from scratches.

Supplied in a nice colour coordinated carrying bag with multilingual box and user manual: GB, FR, IT, ES, DE, PT, GR.