



Product Code	32729
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	130717
EAN/UPC	8023279327298
PARAF	903174959
GTIN	08023279327298

Description

SIRIO SPHYGMO - adult

Classic sphygmomanometer with fabric cuff with velcro closure (green/grey colour).
 Reliable large dial manometer (60 mm) for easy readings, fixed on bulb.
 Scale 0-300 mmHg, graduated every 2 mmHg.
 Supplied in carrying case and box, with multilingual user manual: GB, FR, IT, ES.