



Product Code	32731
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	10049
CND	C9006
EAN/UPC	8023279327311
UMDNS	16156
PARAF	903174961
GTIN	08023279327311
GMDN	16156

Description

BOSTON SPHYGMO - adult - black cuff

Strong durable and elegant palm model.

Aneroid metal zinc-alloy manometer with integral PVC bulb and air release valve for one hand use, provided with a special left and right hand adjustable handle, nylon velcro closure calibrated cuff and 1 tube TPU bladder.

Easy to read: large dial Ø 60 mm.

Unique push button air release valve offers precise mini-deflation.

Special PVC patented ribber bulb.

Latex free.

Supplied with a nice nylon carrying bag.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, Arabic.