

R1 SHOCK PROOF SPHYGMOMANOMETER with one-piece cuff - 1250-150



Product Code	32732
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1039287
CND	C9006
EAN/UPC	4045396178441

Description

R1 SHOCK PROOF SPHYGMO with one piece cuff - 1250-150

Exact precision that lasts a lifetime.

- shock-proof up to a falling height of 120 cm
- valve base and air duct with metal precision valve
- two-component ABS casing
- metal threaded tube connection at the top of the casing for fast cuff exchange and ergonomical operation
- maximum error tolerance of ± 3 mmHg
- supplied in zipper bag

Non-resistant opening of air release valve ensures precise and comfortable operation

Scale: \varnothing 52 mm black with red fluorescent pointer

Fast, ergonomic inflation due to optimized bulb design with integrated spoon

Weight: 173 g

Made in Germany.