



Product Code	32733
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

EAN/UPC 6970384254563

Description

DALLAS SPHYGMOMANOMETER - trolley model

Provided with rear plastic basket and trolley adjustable in height, with 5 plastic castors Ø 5 cm.

All DALLAS sphygmometers are provided with extra large dial plate (Ø 147 mm) for easy reading, black nylon cuff with D-ring and 2 tube adult latex bladder.

Scale 0-300 mmHg.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

Multilingual box (32749 only): GB, FR, IT, ES, PT, DE, GR, Arabic.

Rear plastic basket for cuff storage.

Spiral PVC tube (1.8 m)

Dual-use air-release valve

Adjustable height 95 to 135 cm

5 castors with brake