



Product Code	32736
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2944139
CND	Z1203020302
EAN/UPC	8023279327366
UMDNS	16156

Description

TOKYO 2 SPHYGMO - adult - black cuff

Professional sphygmomanometer with large Ø 65 mm dial and button type air release valve. ABS light plastic case and extra shock-proof ring for long durability. Supplied with black nylon calibrated cuff (54x14.5 cm) with D-ring, black 1-tube PVC, latex free bladder and black, latex free PVC bulb in a nylon carrying bag and cardboard box.

Multilingual box and manual: GB, FR, IT, ES, DE, PT, GR, Arabic.

Flat bottom bulb with three handling modes for right and left users