

ROMA SPHYGMOMANOMETER SET - multi-cuff



Product Code	32738
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	18161
CND	C9006
EAN/UPC	8023279327380
UMDNS	16156
PARAF	903174858
GTIN	08023279327380
GMDN	16156

Description

ROMA KIT - Sphygmo with 3 cuffs: one pediatric (35x11 cm), one adult (53x14.5 cm) and one adult large (62x17.5 cm).

Classic, lightweight, reliable palm model aneroid sphygmo with certified manometer ABS housing, integral PVC bulb and air release valve for one-hand use, provided with handle, blue nylon, velcro closure, calibrated cuff and 1 tube TPU bladder.

Easy to read: large dial Ø 60 mm, with rubber protective ring.

Supplied in a nice nylon carrying bag with external velcro pocket and cardboard box.

Made with no natural rubber or latex.