

# "G7" HEINE SPHYGMOMANOMETER



Product Code	32739
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
--------	-----

RDM (NSIS)	85576
CND	C9006
EAN/UPC	4053755136346
UMDNS	16-156
PARAF	987475213
GMDN	16156

## Description

M-000.09.232 HEINE GAMMA G7 - adult - Adult cuff (29-41 cm) - Zipper case

HEINE GAMMA® G7 Exclusive air-release valve provides fine and precise air release rate adjustment, quick deflation, and extreme durability. Upper housing made of aluminium alloy with a shock absorbing ring protecting indicator and glass. Ergonomic design for comfort and ease of use.

The HEINE GAMMA series is designed for life-long accuracy, every model is shock-resistant, latex-free and features a pressure-proof manometer protected by a micro-filter.

- high quality indicator, precisely calibrated, durable and shock resistant
- large Ø 56 mm black-on-white scale for easy, accurate readings
- housing is moulded of shockproof plastic with shock absorbing, anti-slip coating
- large insufflation bulb (50% larger) for fast cuff inflation
- adjustable spoon for both right and left-handed operation
- high-quality, washable single-tube cuff with velcro fastener

Made in Germany.