

# BIG BEN SPHYGMOMANOMETER - trolley - adult one piece cuff



|               |                |
|---------------|----------------|
| Product Code  | 32740          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I M            |

|         |               |
|---------|---------------|
| EAN/UPC | 4045396188853 |
|---------|---------------|

## Description

RIESTER BIG BEN® - trolley Ø 54 cm - adult one piece cuff

Aneroid sphygmomanometer, made of resistant plastic, with large dial (147 mm) and white pointer on black dial plate for easy-to-read measurements from 0 mmHg to 300 mmHg.  
Basket on back for cuff storage after use.  
Velcro-closure black cuff 2 tubes and long coil tubing. Latex free.

Adjustable height 90-140 cm  
New trolley Ø 54 cm

Made in Germany