

"G5" HEINE SPHYGMOMANOMETER without carrying bag



Product Code	32741
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
RDM (NSIS)	85576
CND	C9006
EAN/UPC	4053755136667
UMDNS	16-156
PARAF	987475201
GMDN	16156

Description

M-000.09.561 HEINE GAMMA G5 - adult

HEINE GAMMA® G5

Exclusive air-release valve provides fine and precise air release rate adjustment, quick deflation, and extreme durability.

The HEINE GAMMA series is designed for life-long accuracy, every model is shock-resistant, latex-free and features a pressure-proof manometer protected by a micro-filter.

high quality indicator, precisely calibrated, durable and shock resistant

large Ø 56 mm black-on-white scale for easy, accurate readings

housing is moulded of shockproof plastic with shock absorbing, anti-slip coating

large insufflation bulb (50% larger) for fast cuff inflation

adjustable spoon for both right and left-handed operation

high-quality, washable single-tube cuff with velcro fastener

Made in Germany.