



Product Code	32743
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	94919
CND	C9006
EAN/UPC	4045396102675

Description

RIESTER BIG BEN® - wall - Rotation 130°, 65° right, 65° left

Aneroid sphygmomanometer, made of resistant plastic, with large dial (147 mm) and white pointer on black dial plate for easy-to-read measurements from 0 mmHg to 300 mmHg.

Basket on back for cuff storage after use.

Velcro-closure black cuff 2 tubes and long coil tubing. Latex free.

Made in Germany