

# DENVER SPHYGMOMANOMETER - trolley - nylon cuff, TPU bladder



Product Code	32744
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	3016168
EAN/UPC	8023279327441
UMDNS	16156
PARAF	987471760

## Description

### DENVER SPHYGMOMANOMETER - trolley - adult cuff

Trolley sphygmomanometer, light but with high stability and 360° swivelling castors.

Adjustable height (90 to 135 cm) with metal basket.

Square-shaped dial 145x145 mm.

Nylon cuff, 2-tube adult TPU bladder, PVC bulb, dual use air-release valve with spring, spiral PVC tube (1.8 m).

Multilingual manual: GB, FR, IT, ES, PT, DE, GR.