



| | |
|---------------|----------------|
| Product Code | 32745 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I M |

| | |
|------------|----------------|
| UK-REP | Yes |
| CH-REP | No |
| RDM (NSIS) | 10049 |
| CND | C9006 |
| EAN/UPC | 8023279327458 |
| UMDNS | 16156 |
| GTIN | 08023279327458 |
| GMDN | 16156 |

Description

DAYTON SPHYGMO - desk/wall model

Provided with big dial (Ø 113 mm) for easy reading.
Basket with handle easy to move and carry.
Nylon cuff with D ring. Wall fixing.
Size: 130x130x135 mm - Weight: 180g
Measuring range: 0-300 mmHg ± 3 mmHg.

Large basket with handle for cuff

PVC bladder and bulb

LATEX FREE