



Product Code	32745
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	No
RDM (NSIS)	10049
CND	C9006
EAN/UPC	8023279327458
UMDNS	16156
GMDN	16156

## Description

DAYTON SPHYGMO - desk/wall model

Provided with big dial (Ø 113 mm) for easy reading.  
Basket with handle easy to move and carry.  
Nylon cuff with D ring. Wall fixing.  
Size: 130x130x135 mm - Weight: 180g  
Measuring range: 0-300 mmHg ± 3 mmHg.

Large basket with handle for cuff

PVC bladder and bulb

LATEX FREE