



Product Code	32746
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1151937
CND	C9006
EAN/UPC	8023279327465
PARAF	903175065

Description

RIESTER BIG BEN® - trolley Ø 54 cm - velcro adult cuff

Aneroid sphygmomanometer, made of resistant plastic, with large dial (147 mm) and white pointer on black dial plate for easy-to-read measurements from 0 mmHg to 300 mmHg.

Basket on back for cuff storage after use.

Velcro-closure black cuff 2 tubes and long coil tubing. Latex free.

Adjustable height 90-140 cm

New trolley Ø 54 cm

Made in Germany.