



Product Code	32747
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	Yes

RDM (NSIS)	10049
CND	C9006
EAN/UPC	8023279327472
PARAF	987471758
GTIN	08023279327472

## Description

DALLAS SPHYGMOMANOMETER - trolley model

Provided with rear plastic basket and stainless steel trolley adjustable in height, with 5 plastic castors Ø 5 cm.  
All DALLAS sphygmometers are provided with extra large dial plate (Ø 150 mm) for easy reading, black nylon calibrated cuff and latex quality inflation system and bladder.  
Scale 0-300 mmHg.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

5 castors (1 with brake)  
Stainless steel trolley  
Base Ø 60 cm