



Product Code	32749
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	10049
CND	C9006
EAN/UPC	8023279327496
UMDNS	16156
PARAF	903627735
GTIN	08023279327496
GMDN	16156

Description

DALLAS SPHYGMOMANOMETER - wall model

Provided with rear plastic basket for cuff storage.

All DALLAS sphygmomanometers are provided with extra large dial plate (Ø 147 mm) for easy reading, black nylon cuff with D-ring and 2 tube adult latex bladder.

Scale 0-300 mmHg.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

Multilingual box (32749 only): GB, FR, IT, ES, PT, DE, GR, Arabic.