



|               |                |
|---------------|----------------|
| Product Code  | 32750          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I M            |
|               |                |
| UK-REP        | Yes            |
| CH-REP        | Yes            |
|               |                |
| EAN/UPC       | 4045396001886  |
| PARAF         | 987475074      |

Description

Exacta ® SPHYGMO

High quality precision sphygmomanometer twin-tube model with black lacquered metal casing

- life-long precision thanks to practically non-ageing and pressure-resistant measuring system up to 600 mm Hg
- stainless metal clamp for hooking the pressure gage to the cuff
- precision air release valve with fine adjustment and wear-free
- microfilter protects air release valve and measuring system
- linear scale Ø 49 mm from 0 to 300 mmHg
- latex bulb for pumping up the cuff
- maximum error tolerance ±3 mmHg

Supplied in vinyl zipper case.

Made in Germany.