



Product Code	32750
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396001886
PARAF	987475074

Description

Exacta ® SPHYGMO

High quality precision sphygmomanometer twin-tube model with black lacquered metal casing

- life-long precision thanks to practically non-ageing and pressure-resistant measuring system up to 600 mm Hg
- stainless metal clamp for hooking the pressure gage to the cuff
- precision air release valve with fine adjustment and wear-free
- microfilter protects air release valve and measuring system
- linear scale Ø 49 mm from 0 to 300 mmHg
- latex bulb for pumping up the cuff
- maximum error tolerance ± 3 mmHg

Supplied in vinyl zipper case.

Made in Germany.