



Product Code	32751
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396178731
PARAF	987474931

## Description

E-MEGA® SPHYGMO - black - adult

Disinfectable, washable and latex free one-piece cuff.

- fast inflation of the cuff due to ergonomic spoon design
- lifelong precision due to almost non-aging precision movement, with specially tempered copper-beryllium diaphragm for pressure-loading capacity up to 600 mmHg
- maximum error tolerance  $\pm 3$  mmHg
- comfort - air release valve, wearfree, with fine adjustment

Made in Germany.

High-polished manometer surface, easy to disinfect  
Large dial  $\varnothing$  63 mm  
Supplied in zipper bag