

PRECISA N SPHYGMO WITH STETHOSCOPE - plastic/metal - 1 tube



Product Code	32754
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	118366
CND	C9006
EAN/UPC	4045396112995
PARAF	950920708

Description

PRECISA® N with stethoscope - plastic/metal - 1 tube

High quality, reliable sphygmomanometers with precision that lasts a lifetime.

- microfilter protects air release valve and measuring system
- specially tempered copper-beryllium diaphragm, with precision movement and almost non-aging
- diaphragm pressure-loading capacity up to 600 mmHg
- metal valve base
- ergonomic spoon-shaped grip made of stainless steel
- linear plastic/metal, scale from 0 to 300 mmHg (easily legible)

Supplied with calibrated nylon cuff with velcro in synthetic leather bag.

Made in Germany.