



Product Code	32756
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396110786
PARAF	987475643

Description

PRECISA® N - aluminium - 2 tubes

High quality, reliable sphygmomanometers with precision that lasts a lifetime.

- microfilter protects air release valve and measuring system
- specially tempered copper-beryllium diaphragm, with precision movement and almost non-aging
- diaphragm pressure-loading capacity up to 600 mm Hg
- metal valve base
- ergonomic spoon-shaped grip made of stainless steel
- linear, tilted aluminum scale, from 0 to 300 mm Hg (easily legible)

Supplied with calibrated nylon cuff with velcro in synthetic leather bag.

Made in Germany.