

PRECISA N SHOCK PROOF SPHYGMO - aluminium - 1 tube - adult



Product Code	32757
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2812201
CND	C9006
EAN/UPC	4045396211681
PARAF	987475821

Description

PRECISA® N SHOCK PROOF SPHYGMO

Shock-resistant aneroid sphygmomanometer with same features as Precisa® N.

Shock-resistant from a height of 1 m, according to DIN EN 81060-1:2012-08, to prevent impact or falls on the floor from compromising measurement accuracy.

High quality, reliable sphygmomanometers with precision that lasts a lifetime.

- microfilter protects air release valve and measuring system
- specially tempered copper-beryllium diaphragm, with precision movement and almost non-aging
- diaphragm pressure-loading capacity up to 600 mm Hg
- metal valve base
- ergonomic spoon-shaped grip made of stainless steel

Supplied with calibrated nylon cuff with velcro in synthetic leather bag.

Linear aluminium scale Ø 63 mm.

Made in Germany.