

RIESTER SPHYGMOTENSIOPHONE - 1380



Product Code	32760
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M
EAN/UPC	4045396001930
PARAF	903175356

Description

SPHYGMOTENSIOPHONE - 1380

Extreme precise two-tube sphygmomanometer with refined functionality.
High quality cotton adult cuff finished in artificial leather is provided with a metal clip through which to pull the cuff.

Supplied in zippered nylon bag.

- precision air release valve with fine adjustment and wear-free
- microfilter protects air release valve and measuring system
- tempered copper-beryllium diaphragm almost non-aging
- diaphragm pressure-loading capacity up to 600g mmHg
- maximum error tolerance of ± 3 mmHg
- easy to read scale $\varnothing$ 49 mm from 0 to 300 mmHg

Made in Germany.