

# PRECISA N SHOCK PROOF SPHYGMO - aluminium - 1 tube - pediatric - 1364-129



|               |                |
|---------------|----------------|
| Product Code  | 32761          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I M            |

|            |               |
|------------|---------------|
| UK-REP     | Yes           |
| CH-REP     | Yes           |
| RDM (NSIS) | 2934931       |
| CND        | C9006         |
| EAN/UPC    | 4045396211704 |

## Description

### PRECISA® N SHOCK PROOF SPHYGMO - pediatric - 1364-129

Shock-resistant aneroid sphygmomanometer with same features as Precisa® N.

Shock-resistant from a height of 1 m, according to DIN EN 81060-1:2012-08, to prevent impact or falls on the floor from compromising measurement accuracy.

High quality, reliable sphygmomanometers with precision that lasts a lifetime.

- microfilter protects air release valve and measuring system
- specially tempered copper-beryllium diaphragm, with precision movement and almost non-aging
- diaphragm pressure-loading capacity up to 600 mm Hg
- metal valve base
- ergonomic spoon-shaped grip made of stainless steel

Supplied with calibrated nylon cuff with velcro in synthetic leather bag.

Linear aluminium scale Ø 63 mm.

Made in Germany.