

EMERGENCY MULTICUFF SPHYGMO KIT with 5 cuffs



Product Code	32771
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1711341
CND	C9006
EAN/UPC	6945799153120
PARAF	987472622

Description

EMERGENCY SPHYGMO KIT - 5 cuffs

Palm type sphygmomanometer with orange silicone shockproof ring and orange cup of air-release valve. PVC tube and 1 tube TPU bladder. Supplied in a plastic carrying case with 5 cuffs: XS infant (13-19 cm), S child (18.4-26.7 cm), M adult (25.4-40.6 cm), L adult large (33.3-51 cm), XL thigh (40.6-66 cm).

Multilingual manual and box: GB, FR, IT, ES.