



Product Code	32775
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1334503
CND	Z1203020599
EAN/UPC	5060006581311

Description

UM-102A MERCURY FREE SPHYGMOMANOMETER - desk

AND professional highly advanced digital blood pressure monitor

- auscultatory blood pressure measurement
- digital readout for pulse rate
- cuff holder and grip for carrying
- durable and chemical resistant body
- full ABS housing lightweight body
- adjustable display column for easier analysis of results
- anti-glare display minimizes reflected light

Supplied with adult cuff 22-32 cm

Made in Japan.

Technical Specifications

Size: 360x96x66 mm

Measuring range

- pressure: 0-300 mmHg \pm 3 mmHg

- pulse: 40-180 beats/minute \pm 5%

Power supply: 2x1.5V AA alkaline batteries