



Product Code	32777
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1671668
CND	Z1203020501
EAN/UPC	8023279327779
UMDNS	16-157
PARAF	982172506

Description

JOLLY BLOOD PRESSURE MONITOR

Professional Blood Pressure monitor for medical use.

- large easy-to-read illuminated LCD display 86x74 mm
- WHO blood pressure indicator
- arrhythmia detection (IHB)
- 4 user memories, each with 30 measured values, time and date
- average of all saved measured values and morning and evening blood pressure of the last 7 days
- time and date
- automatic switch-off
- Medical Devices Class IIa

Supplied with adult cuff, soft storage case, batteries and multilingual manual GB, FR, IT, ES, PT, DE, PL, GR.

Technical Specifications

Measurement system: Oscillometric method
 Measurement range: Pressure: 0300 mmHg \pm 3 mmHg
 Pulse: 40180 pulse/min. \pm 5%
 Power supply: 4x1.5 V AAA batteries
 Size: 150x95xh 41 mm
 Weight (without batteries): 69.5 g approx 235 g