



Product Code	32778
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2656013
CND	C9006
EAN/UPC	6970384250107
PARAF	950079956

Description

DOMINO DIGITAL BLOOD PRESSURE MONITOR

Professional Blood Pressure monitor.

- extra large LCD display with innovative design
- date/time/ambient temperature (°C/°F) display
- IHB Irregular heartbeat detector
- WHO blood pressure classification indicator
- 120 memories (60 sets for 2 users)
- alarm clock reminding function

Supplied with adult cuff, storage case, 4xAA 1.5V batteries and 1xCR2025 button battery.

Multilingual manual and box: GB, FR, IT, ES, PT, DE, NL, SK, RO, HR, GR, BG

Extra large display 115x79.3 mm

Technical Specifications

Measuring range: Pressure - 40~260 mmHg / Pulse: 40~160 beats/min.

Accuracy: Pressure ± 3 mmHg Pulse: $\pm 5\%$

Power supply: 4x "AA" batteries or AC adaptor (option)

Size - weight: 158x120xh 127 mm - 490 g