



Product Code	32786
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	237582
EAN/UPC	8023279327861
GTIN	08023279327861

## Description

### GIMA AUTOMATIC WRIST B.P. MONITOR

Reliable digital blood pressure monitor with following features:

- one touch automatic inflation and deflation
- large LCD display with date and time
- 2 person use / 60 memories each
- power saving system

Multilanguage box and manual (GB, IT, FR, ES, PT, GR, DE, Arabic)

## Technical Specifications

Measuring method: oscillometric  
 Measuring range:  
 Pressure 20~280 mmHg  
 Pulse 40~195

Accuracy:  
 Pressure:  $\pm 3$  mmHg  
 Pulse:  $\pm 5\%$  of reading

Power supply: 2x "AAA" batteries  
 Dimensions: 80x60x75 mm  
 Weight: approx 130 g (including cuff, ex batteries)  
 Wrist circumference: 13.5 - 19.5 cm