

# LED DISPLAY MERCURY-FREE SPHYGMOMANOMETER - desk



Product Code	32797
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2572524
CND	C900399
EAN/UPC	8023279327977

## Description

LED DISPLAY MERCURY-FREE SPHYGMOMANOMETER - desk - light blue

LED mercuryfree sphygmomanometers provide reliable and precise results.

- manual pressurization by bulb
- manual depressurization by air release valve
- LED column and display

Column flashes when pressure is more than 315 mmHg.

Suitable for adult and pediatric use.

Supplied with PVC bladder and bulb, nylon adult cuff (25-38 cm) with D-ring, user manual in GB, FR, IT, ES, PT, DE, PL, CZ, SK, HU, RO, GR, HR, SI, EE, Arabic.

## Technical Specifications

Pressure range: 0-300 mmHg,  $\pm 3$  mmHg

Pulse rate range: 30-200bpm,  $\pm 5\%$

Power supply: 3xAA batteries, DC 4.5V

Weight: 584 g without battery and cuff

Protection: Type BF applied part; IPx21