



Product Code	32804
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	237582
CND	Z1203020501
EAN/UPC	8023279328042
UMDNS	16157
PARAF	987471822

## Description

Professional desk Blood Pressure monitor for hospital and medical use.

- Extra large LCD display with innovative design
- Date/time/ambient temperature (°C/°F) display
- IHB Irregular heartbeat detector
- 120 memories (60 sets for 2 users)
- Alarm clock reminding function
- WHO blood pressure classification indicator

## Technical Specifications

- Measuring range:
  - Pressure: - 40 ~ 260 mmHg
  - Pulse: 40 ~ 160 bpm
- Accuracy:
  - Pressure  $\pm$  3 mmHg
  - Pulse:  $\pm$  5%
- Power supply: 4x "AA" batteries or AC adaptor (option)
- External dimension: 158 x 120 x 127 mm
- Weight: 490 g