

PREMIUM NEWBORN CUFF + 1 TUBE TPU BLADDER



Product Code	32817
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2563718
CND	C99
EAN/UPC	8023279328172
GTIN	08023279328172

Description

- PREMIUM CUFFS dark grey - nylon - TPU bladder - velcro fastening
 Made of high quality nylon 210D.
 Durability test up to 10,000 cycles.
- One tube cuff for Sirio, London, Boston, Roma 1T, Tokyo, Madrid, Ri-san, Precisa, Minimus...
 - Model: Newborn
 - Size cm: 5.5x18.5
 - Arm circumference: 9.5 - 13.5 cm