

PREMIUM ADULT CUFF + 1 TUBE BLADDER - dark grey



Product Code	32830
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512519
CND	C99
EAN/UPC	8023279328301
PARAF	903170520
GTIN	08023279328301
GMDN	34978

Description

PREMIUM CUFFS dark grey - nylon - TPU bladder - velcro fastening
Made of high quality nylon 210D.
Durability test up to 10,000 cycles.

- One tube cuff for Sirio, London, Boston, Roma 1T, Tokyo, Madrid, Ri-san, Precisa, Minimus...
- Model: Adult
- Size cm: 54x15
- Arm circumference: 25.4 - 40.6 cm