

PAEDIATRIC CUFF + 2 TUBES BLADDER 35x11 cm - nylon



| | |
|---------------|----------------|
| Product Code | 32847 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2830555 |
| CND | C99 |
| EAN/UPC | 8023279328479 |
| PARAF | 903170672 |
| GMDN | 34978 |

Description

Calibrated - latex free - nylon cuffs with TPU bladder and PVC tube.

- Two tubes cuff for Minor, Yton line, Roma 2T, Dallas, Big Ben, Empire, Dayton...
- Colour: Blue
- Model: Pediatric
- Size cm: 35x11
- Arm circumference: 19 - 26 cm