

PAEDIATRIC CUFF + 2 TUBES BLADDER 35x11 cm - nylon



Product Code	32847
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2830555
CND	C99
EAN/UPC	8023279328479
PARAF	903170672
GTIN	08023279328479
GMDN	34978

Description

Calibrated - latex free - nylon cuffs with TPU bladder and PVC tube.

- Two tubes cuff for Minor, Yton line, Roma 2T, Dallas, Big Ben, Empire, Dayton...
- Colour: Blue
- Model: Pediatric
- Size cm: 35x11
- Arm circumference: 19 - 26 cm