

ADULT CUFF + 2 TUBES BLADDER 49x15 cm - nylon black



Product Code	32850
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512519
CND	C99
EAN/UPC	8023279328509
GTIN	08023279328509
GMDN	34978

Description

Calibrated - latex free - nylon cuffs with TPU bladder and PVC tube.

- Two tubes cuff for Minor, Yton line, Roma 2T, Dallas, Big Ben, Empire, Dayton...
- Colour: Black
- Model: Adult
- Size cm: 49x15
- Arm circumference: 26 - 40 cm