


PRESS-SAN HYGIENIC PROTECTION OF SPHYGMOMANOMETER



Product Code	32858
--------------	-------

Unit of sale	box of 10 pcs
--------------	------------------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279328585
---------	---------------

PARAF	974039974
-------	-----------

Description

PRESSAN® HYGIENIC PROTECTION OF SPHYGMO - box of 10

A practical hygienic accessory for the sphygmomanometer to protect patient's skin whilst measuring the blood pressure.

Made of a white non-woven fabric.

Ensures patient's comfort and hygiene.

Can fit any kind of sphygmomanometer or pressure holter.

Made in Italy.