

MICROLIFE WATCH BP HOME A B.P.MONITOR + SOFTWARE



Product Code	32868
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A
UK-REP	No
CH-REP	No
RDM (NSIS)	2527514
CND	Z1203020501
EAN/UPC	4719003310509

Description

MICROLIFE Watch BP® HOME A + SOFTWARE (Professional)

- early detection of Atrial Fibrillation for stroke risk prevention
- assists patients in following procedures for accurate home blood pressure measurement between doctor visits
- usual mode: 3 reading average mode with 250 data memories
- diagnostic mode: 7-day self-measurement program that follows ESH/AHA/BIHS guidelines for accurate measurements at home, providing the doctor with trustworthy results
- user-friendly software to transfer data to PC via USB
- 3 consecutive measurements for high accuracy
- clinically validated also for children (over 12 years), diabetics, and End-Stage Renal Disease (ESRD) patients

Supplied with soft cuff 22-32 cm (M size), software (GB), quick start guide, 4x1.5V AA batteries, carrying bag, and manual (GB, IT, FR, DE).

Other languages available on gimaitaly.com (ES, DK, FI, SE, NL, NO, GR, UA, RU).

Size: 150x100x50 mm

Weight: 385 g

AFIB detection
5 years warranty